Flash Sterilization Trouble Shooting Guide

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Trouble Shooting Guide

The purpose of this trouble shooting guide is to help you identify the possible cause(s) of sterilization processing failures in 270°F gravity displacement (flash) sterilizers. It is designed to highlight potential problem areas through the use of 3 worksheets and the explanations provided here. By focusing on these problem areas, you can reduce the amount of time it takes to put your machine back into service.

Section I Pinpointing the Problem

(Use Worksheet I)

Worksheet I is designed to help you review your records to determine if the problem is widespread.

It may be helpful to fill in the information on Worksheet I before reading further.


As you review the information on Worksheet I, look for situations where patterns appear. For instance, if you have two machines on the same steam line that are producing positive biological results, steam delivery may be a problem. You may want to involve your engineering department. On the other hand, if the problem only occurs when a specific operator is using the machines, you may wish to investigate how that person operates the sterilizer. When you have a problem with only one machine, and you can focus on that machine to find a solution, then go directly to Worksheet II.

Section II Profiling a Problem Machine

(Use Worksheet II)

A normal emergency (flash) gravity displacement autoclave cycle at 270°F should operate like this:

Figure 1

Time / Temperature

Come-up is the time it takes for the chamber to reach the temperature of the incoming saturated steam. The exposure time of the cycle should not start until all air has been removed and the temperature in the drain line stops rising.

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
Pressure

Saturated steam is assured when the temperature and pressure have the following relationship:

Figure 2

Temperature (°F)	Pressure (p.s.i.)
268.2	26
269.7	27
271.2	28
272.6	29
273.9	30
275.3	31
276.6	32
278.0	33
279.3	34
280.5	35

This relationship should be maintained throughout the exposure time.

Cycle Analysis

Worksheet II is intended to focus on the exact cycle conditions for each problem machine. You will most likely find that each machine has its own unique cycle pattern; however, each machine may vary from cycle to cycle.

In order to detect cycle variability, it will be necessary to run several cycles in each problem machine and record the results on Worksheet II.

As you review your information on Worksheet II, you may observe that your *timer* is activated before the temperature has reached its *highest point*. If this happens, it indicates that the time of the exposure portion of the cycle started prematurely. Usually *Come-up times* are usually between 1 and 1.5 minutes.

Machines with come-up times of less than 1 minute usually show the timing to be started before the temperature has reached its highest point. Remember, a mixture of air and steam will not be as hot as saturated steam. Therefore, if the temperature in your chamber continues to rise when your timer is running, your *exposure time* has started prematurely and the chamber does not contain saturated steam.

If your temperature (as recorded by various sensors) and pressure (recorded by the chamber pressure gauge) do not agree with the values listed in Figure 2 when the highest temperature is recorded, this can indicate either:

- 1) gauge error
- 2) steam quality variation, or
- 3) a combination of gauge error and steam quality variation

Variation may also show up when you compare the cycle pattern of a machine when it is full vs. the cycle pattern when it is empty. These variations are to be expected and may account for intermittent positive biological results.

There are a number of ways the cycle variation can be managed as outlined in Section IV.

Section III Reviewing Spore Testing Procedure (Use Worksheet III)

The biological monitor selected to challenge these sterilizers should be of consistently high resistance so that the indicator can detect variations in the sterilizers' performance.

It may be helpful to fill in the information on Worksheet III before reading further.

Just as it is important to assure the consistency of the sterilizer cycle conditions, it is important to be consistent in the use of your biological indicator.

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Testing at the same time every day is a common practice. However, if testing is done "first thing in the morning", some sterilizers may have condensed steam in the steam line from standing unused all night. A "blank" cycle should be run on each machine if a "first thing in the morning" schedule is used to clear steam condensate from the line.

The test cycle should have the same conditions as the cycles you use for routine reprocessing of items. In fact, to maximize the usable information given by your biological indicator, you should include a biological indicator with every load being sterilized, or with a test load that is representative of your actual loads.

All biological indicators should be used according to the recommendations of the manufacturer. A review of manufacturer's instructions may be helpful.

If all information recorded on the worksheets does not readily indicate a problem but your problems persist, you may wish to have that specific lot of indicators evaluated by their manufacturer.

Section IV Solving the Problem

Finding solutions to sterilization problems is a complex process. There are many variables to consider. Using the worksheets will help to bring these variables into focus. There are a number of resources available to help you in this investigation: hospital maintenance personnel, sterilizer repairmen, biological indicator manufacturers, and professional organizations. Each resource will look at the problem from a different perspective, but it will be up to you to ask the right questions and to make the decisions.

Usually there is no single right answer. You have to weigh the alternative solutions against their respective drawbacks to determine the one that is most appropriate for your situation. For instance, you may

discover that you have one or more sterilizers that vary greatly from cycle to cycle. Your options may be: 1) to buy new machines; 2) to have your machines serviced once a week; or 3) to extend the cycle time to overcome the machines' variability.

Each of the above options could solve the problem of intermittent positive biological results. Buying new machines, or having them serviced each week could prove to be very expensive. A spore test with a shorter kill time may reduce your ability to detect potentially serious sterilization failures early. Extending cycle times would require you to change departmental procedures and retrain the staff.

The process of solving sterilization process failures is complex and the choices are difficult. Careful investigation of the possible solutions is the most effective way to make the cost-effective, ethical and professional decisions.

Section V Maintenance: Problem Prevention

Preventive maintenance on equipment is one of the most effective means of reducing sterilization failures. The list of activities below are very important to the overall consistency and reliability of sterilizer operation:

- cleaning drain lines
- calibrating all gauges
- checking door seals
- verifying that the pressure in the steam lines when all machines are on is the same for all machines
- checking for steam trap maintenance.

Your hospital maintenance personnel and sterilizer repairmen will be able to help you plan and evaluate your maintenance program.

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3M Canada Inc.